





Betaflex Titan



NI-FREE*

Order No.	Profil profile	ø mm	ø inch	Inhalt contents	Beschreibung description
CE 0297 220-4055		0,41 x 0,56	.016 x .022	10	Betaflex Titan Drähte - gerade à 380 mm = 15" Länge
220-4563		0,46 x 0,64	.018 x .025	10	Betaflex titanium wires - straight of 380 mm = 15" length
220-4081		0,8	.032	10	Betaflex Titan Drähte - gerade à 380 mm = 15" Länge, für Pendulum Appliance Betaflex titanium wires - straight of 380 mm = 15" length, for Pendulum Appliance




Ni-Free Products

Bukkalröhrchen zum Kleben – Nickelfrei Buccal tubes for direct bonding – Ni-Free

CE 0297



Pad sandgestrahlt
Sandblasted pad

OK / Maxillary	Zahn tooth	Torque	Angulation	distaler Offset distal offset	auf/on Pad (10 Stück pro Packung/pcs. per pack)		
					rechts/right	links/left	
occlusal headgear tube  convertible	$\frac{6}{ } \frac{6}{ }$	-14°	0°	14°	.018" x .025"	741M0704	740M0704
					.022" x .028"	741M0702	740M0702
	$\frac{7}{ } \frac{7}{ }$	-14°	0°	14°	.018" x .025"	741M0804	740M0804
					.022" x .028"	741M0802	740M0802
UK / Mandibular	$\frac{6}{ } \frac{6}{ }$	-30°	0°	4°	.018" x .025"	741M1704	740M1704
					.022" x .028"	741M1702	740M1702
	$\frac{7}{ } \frac{7}{ }$	-30°	0°	4°	.018" x .025"	741M1804	740M1804
					.022" x .028"	741M1802	740M1802

100T2000


CE



7



0,8



10

Universal Dehnschraube wie 100-2000 aus Titan

Universal expansion screw as 100-2000 made of titanium

-	243-2070	0,70	.028	30 m
-	243-2080	0,80	.031	20 m
-	243-2090	0,90	.036	10 m

FORESTANIT® - Ni Free

Laborrollen

rund, nickelfreier Draht*

Laboratory coils

round, nickel free wires*

* < 0,3 % Ni

